Product Catalogue 2011 - 2012

Enzymes and Chemicals for PCR

GeneON GmbH

Hubertusstraße 20
D-67065 Ludwigshafen

Phone: +49-621-5720 864

Fax: +49-621-5724 462

E-Mail: info@geneon.net

Web: www.geneon.net

Content

Enzymes for Molecular Biology

T4 DNA Ligase..... 8- 9

Chemicals for Molecular Biology

Agarose..... 12-13

Proteinase K..... 14-15

IPTG..... 16-17

Bovine Serum Albumin - BSA..... 18-19

UDG (Uracil-DNA Glycosylase)..... 20-21

Enzymes/Chemicals

Enzymes for Molecular Biology 7

Chemicals for Molecular Biology 11

Enzymes for Molecular Biology

- T4 DNA Ligase..... 8-9

T 4 DNA Ligase

Anwendung:

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T 4 DNA Ligase

Applications:

- Cloning of restriction fragments
- joining linkers and adapters to blunt-ended DNA
- gene-(fragments) synthesis.

Description:

T4 DNA Ligase catalyzes the formation of a phosphodiester bonds between 5' phosphate and 3'-hydroxyl-termini in duplex DNA/RNA. This enzyme can join blunt end and cohesive end termini, repairs single stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids.

Concentration: 100 - 200 u/μl

Source:

Isolated from E.coli strain that carries the cloned DNA ligase gene from bacteriophage T4

Unit definition:

One unit is defined as the amount of enzyme required to give 50% ligation of Hind III fragments of lambda DNA in 30 minutes at 16°C at 5' termini concentration of 0.12 μM (300 μg/ml). One Cohesive End Ligation Unit equals 0.015 Weiss units. One Weiss unit equals 67 Cohesive End Ligation Units.

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
Agarose for gel-electrophoresis

Features:

- High separation properties and sharp band patterns
- Easy solubility without foaming
- Excellent optical transparency

Applications:


- DNA: approx. 0.05 kbp – 50 kbp
- RNA: approx. 0.30 kb – 20 kb

Specifications :

- Gelling temperature: $\leq 37\text{ }^{\circ}\text{C}$
- Melting temperature: $\leq 90\text{ }^{\circ}\text{C}$
- Electroendosmosis: ≤ 0.140
- Gel strength (1.5 %): $\geq 2300\text{ g / cm}^2$
- Sulphate content: $\leq 0.10\text{ %}$
- Water content: $\leq 10.0\text{ %}$

Quality Assurance:

- Certified free of DNases and RNases
- No DNA binding
- High lot-to-lot consistency

1% 1x TAE Universal Agarose gel showing separation of λ -DNA digested with *EcoR I/Hind III* (1) and a mixture from pBR328-DNA digested with *Bgl I* and *Hinf I* (2). Data in kbp.

Cat.-no	Description	Amount
604-001	Agarose universal for gel-electrophoresis	100 g
604-005	Agarose universal for gel-electrophoresis	500 g

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Proteinase K

Features:

Proteinase K is a highly active and stable protease with low cutting specificity. The enzyme belongs to the group of subtilisin-related serine proteases and is strongly inhibited by PMSF.

Description:

Proteinase K is a serine protease that exhibits a very broad cleavage specificity. The Protein with a molecular weight 28.900 kDa cleaves peptide bonds adjacent to the carboxylic group of aliphatic and aromatic amino acids. Proteinase K is not inactivated by chelating reagents such as EDTA or detergents such as SDS and is active over a wide range of pH (4-12.5).

Activity: > 30 units/mg protein (hemoglobin, pH 7.5, 37°C)

Unit definition:

One unit is the amount of enzyme which releases at 37°C in 1 min as many folin-positive amino acids and peptides from hemoglobin as 1 µmol of tyrosine.

Usage:

In presence of 0.5 – 1 % SDS Proteinase K inactivates DNases and RNases in eukaryotic and microbiological cell cultures. The use of Proteinase K during lysis of the cells allows the isolation of intact highly-molecular nucleic acids.

Quality:

- purified by chromatography and lyophilized
- RNases: not detectable
- DNases: not detectable
- Exonucleases: not detectable

Cat.-no	Description	Amount
405-002	Proteinase K	200 mg
405-010	Proteinase K	1000 mg

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Applications:

- Blue/white colony screening
- Expression of cloned genes that are under control of the lac promoter

Description:

IPTG (isopropyl-beta-D-thiogalactopyranoside) is a highly stable synthetic analogue of lactose. It inactivates the lac repressor and induces synthesis of beta-galactosidase, an enzyme that promotes lactose utilization. The IPTG is used to induce the expression of cloned genes which are under control of the lac operon. It is used in conjunction with X-Gal to determine the lac phenotype in blue/white colony screening.

Quality:

- Dioxane-free.
- Greater than 99.8 % purity confirmed by HPLC
- Functionally tested in blue/white colony screening

Usage:

The final concentration differs from author to author (0.1 mM (2) - 1 mM (1) IPTG in LB-medium). For preparing a stock solution dissolve 100 mM or 500 mM IPTG in water or buffer, sterilize by filtration through a 0.22 μ m filter and store at -20°C.

Cat.-no	Description	Amount
406-010	IPTG	1 g
406-011	IPTG	5 g
406-012	IPTG	5x5 g

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BSA (Bovine Serum Albumin)

Description:

Bovine Serum Albumin (BSA) is supplied with some restriction enzymes to prevent adhesion of the enzyme to reaction tubes and pipette surfaces. BSA also stabilizes some proteins during incubation. The source of purified BSA is a molecular biology grade, protease free, fraction V starting material. This material is then processed to remove nucleases.

Concentration: 10 mg/ml

Storage Buffer:

- 20 mM KPO₄
- 50 mM NaCl
- 0.1 mM EDTA
- 5% Glycerol
- pH 7.0 (25°C)

Note:

It is do not recommend using this preparation of BSA as a standard for quantitation purposes or size determination.

The BSA should be diluted 1:100 for standard reactions. At this level the buffer from the BSA have a very small effect, only.

Cat.-no	Description	Amount
209-001	Bovine Serum Albumin (BSA)	1,5 ml
209-005	Bovine Serum Albumin (BSA)	4x1,5 ml

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UDG (Uracil-DNA Glycosylase)

Applications:

- site-directed mutagenesis
- as a probe for protein-DNA interaction studies
- Glycosylase mediated single nucleotide polymorphism detection (GMPD)
- SNP genotyping
- Rapid and efficient cloning of PCR products
- Elimination carry-over contamination in PCR

Description:

The Uracil-DNA Glycosylase (UDG, UNG) catalyzes the hydrolysis of the N-glycosylic bond between the uracil and sugar, leaving an apyrimidinic site in uracil-containing single or double-stranded DNA. Shows no activity on RNA. Molecular weight: 25.6 kDa monomer.

Concentration: 20 - 40 u/µl

Unit definition:

One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [3H]-uracil in a 50 µl reaction containing 0,2 µg DNA (104 – 105 cpm/µg) in 30 minutes at 37°C.

Storage Buffer:

20 mM Tris-HCl (pH 8.0); 50 mM NaCl; 1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol

Reaction Buffer 10X

200 mM Tris-HCl (pH 8.0 at 25°C), 10 mM EDTA, 10 mM DTT. Incubation of 1X buffer at 37°C.

Note:

At temperatures below 55 °C, the activity of uracil-DNA glycosylase is partially restored.

Cat.-no	Description	Amount
111-005	Uracil-DNA Glycosylase (UDG)	200 units
111-025	Uracil-DNA Glycosylase (UDG)	1000 units

● GeneON ●