# Product Catalogue 2011 - 2012

PCR Related

## GeneON GmbH

Hubertusstraße 20  
D-67065 Ludwigshafen

Phone: +49-621-5720 864

Fax: +49-621-5724 462

E-Mail: [info@geneon.net](mailto:info@geneon.net)

Web: [www.geneon.net](http://www.geneon.net)

## Content

### PCR related products

Maximo Taq DNA Polymerase .....	8-9
Maximo Taq DNA Polymerase (2X pre-mix, ready-to-use) .....	10-11
Maximo BLUE-Taq DNA Polymerase (ready-to-load) .....	12-13
Maximo DFS-Plus Taq DNA Polymerase DNA-free .....	14-15
Maximo M-SuperhotTaq DNA Polymerase (qPCR) .....	16-17
m-AntiTaq .....	18-19
Maximo P-SuperhotTaq DNA Polymerase (qPCR) .....	20-21
Pfu/Psp DNA Polymerase .....	22-23
Pfu/Psp RED-DNA Polymerase (ready-to-load) .....	24-25
Pfu/Psp DNA Polymerase (2X Pre-mix, ready to use) .....	26-27
Maximo Tth Polymerase .....	28-29
UDG (Uracil-DNA Glycosylase) .....	30-31

### PCR-Mastermixes for Real-time PCR (qPCR)

#### Colorless-Line

SuperHot qPCR Master Mix RT .....	34-35
SuperHot qPCR Master Mix HY .....	36-37
SuperHot qPCR Master Mix HS .....	38-39

#### Safe-Line

##### with EvaGreen

Eva-Master Mix E1 .....	40-41
Eva-Master Mix E2 .....	42-43
Eva-Master Mix E3 .....	44-45
Eva Master Mix E4 .....	46-47

##### without stain DLP (for dual label probes)

Master Mix DLP1 .....	48-49
Master Mix DLP2 .....	50-51
Master Mix DLP3 .....	52-53
Master Mix DLP4 .....	54-55

#### Lyo-Line

L-Mastermix for lyophilisation (low glycerol content) .....	56-57
---	-------

<b>LightCycler Capillaries .....</b>	<b>58-61</b>
--------------------------------------	--------------

## Content

### PCR Mastermixes for standard PCR

Maximo Taq DNA Polymerase (2X pre-mix, ready-to-use) .....	64-65
Maximo Blue-Taq DNA Polymerase (ready-to-load, conc. 1u/μl) .....	66-67
Pfu/Psp RED DNA Polymerase (ready-to-load) .....	68-69
Pfu/Psp DNA Polymerase (2X Pre-mix, ready to use).....	70-71

### Reverse Transcription


AMV Reverse Transcriptase .....	74-75
Ribonuclease Inhibitor / RNase Inhibitor .....	76-77
Reverse (M-MuLV RT) .....	78-79
Maximo Tth Polymerase.....	80-81
Oligo(dT) <sub>15</sub> (1DA) .....	82-83
Random Primers .....	84-85

### Nucleotides Ultrapure

Nucleotide-Mix 40mM (4x10mM each)/Nucleotide Set 100mM (4 x 1 ml).....	88-89
Biotin-11-dUTP (1mM) .....	90-91

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# Maximo Taq DNA Polymerase

## Features:

Maximo Taq DNA Polymerase provides robust PCR performance in a wide range of PCR applications and different templates. Best value in terms of cost per unit.

## Applications:

- Standard / General PCR
- Multiplex PCR
- High-throughput PCR
- Primer extension
- Gene mutation
- T/A cloning

## Description:

Maximo Taq DNA Polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-strand specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94 KDa.

**Concentration:** 5 u/μl

## Unit definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitation material in 30 min at 74 degree.

## Storage Buffer:

25mM Tris-HCl (pH8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, 0.5% Nonident P40, 0.5% Tween 20

## Reaction Buffers supplied with the enzyme:

**10X Buffer:** 500mM KCl, 100mM Tris-HCl, pH 9.0, 1% Triton X-100, 15mM MgCl<sub>2</sub>

Cat.-no	Description	Amount
S 101	Maximo Taq DNA Polymerase	500 units
S 102	Maximo Taq DNA Polymerase	5x500 units
S 103	Maximo Taq DNA Polymerase	20x500 units
S 104	Maximo Taq DNA Polymerase	100x500 units

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## (2X pre-mix, ready-to-use)

### Features:

Maximo Taq DNA Polymerase 2X-preMix provides robust PCR performance in a wide range of PCR applications and different templates. Best value in terms of cost per unit. The optimized mixture of all components reduces pipetting mistakes and ensures repeatable results - every day.

### Applications DFS-Taq Polymerase:

- Standard / General PCR
- optimized for high specificity
- High-throughput PCR, automated pipetting, or plate based PCR
- Gene mutation
- T/A cloning

### Description:

Maximo Taq DNA Polymerase 2X-preMix is optimized and ready-to-use mixture of all components for a successful PCR. Only your primers and your DNA Template has to be added. Maximo Taq DNA Polymerase 2X-preMix contains a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-stranded specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94KD.

**Concentration:** the mixture is 2X concentrated

### Unit definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitation material in 30 min at 74 degree

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(ready-to-load)

## Features:

Maximo Taq-Blue DNA Polymerase provides robust PCR performance in a wide range of PCR applications and different templates. After PCR reaction the enzyme can be loaded to agarose gel, no dye and DNA loading buffer is needed. The enzyme is time- and cost saving, because it includes dye and loading buffer, already.

## Applications:

- Standard / General PCR with visible control
- High-throughput PCR
- Primer extension
- Gene mutation
- T/A cloning

## Description:

Maximo Tag-Blue DNA Polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-stranded specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94 kDa.

**Concentration:** 1 u/μl

## Unit definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitation material in 30min at 74 degree

## Storage Buffer:

25mM Tris-HCl (pH8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, 0.5% Nonident P 40, 0.5% Tween 20, in blue loading buffer

## Reaction Buffers supplied with the enzyme:

**10X Buffer:** 500mM KCl, 100mM Tris-HCl, pH 9.0, 1% Triton X-100, 15mM MgCl<sub>2</sub>  
**MgCl<sub>2</sub>:** 100 mM

Cat.-no	Description	Amount
S 111	Maximo Taq-Blue DNA Polymerase	500 units
S 112	Maximo Taq-Blue DNA Polymerase	5x500 units
S 128	Maximo Taq-Blue DNA Polymerase	20x500 units
S 129	Maximo Taq-Blue DNA Polymerase	100x500 units (or as Bulk)

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## DNA-free

### Features:

DFS-Plus Taq DNA Polymerase provides a new formula in buffers and additives to prevent failures in PCR-applications where inhibitors (e.g. proteins, fat or PS) reduce the performance. The robust enzyme is well suited for sensitive experiments using random primers or bacterial templates. Because of the high sensitivity less than 6 molecules can be detected.

### Application:

Instead of conventionally purified Taq-DNA Polymerase for sensitive PCR reactions, for the detection of bacterial DNA or for applications where inhibitors decrease the performance of regular polymerases

### Reaction conditions:

Same as for conventionally purified Taq-DNA Polymerase.

### Concentration

5 units/µl supplied in 10 mM KPO<sub>4</sub> (pH 7.4 at 25°C), 0.1 mM EDTA, 0.1% Tween 20, 0.1% Triton-X 100 and 50 % (v/v) glycerol.

### Unit Definition

One unit is defined as the amount of enzyme required to incorporate 10 nmoles of dNTP into acid-insoluble material in 30 min at 74°C.

### Reaction Buffers provided:

Ammonium-Reaction buffer (10X) "incomplete"  
Ammonium-Reaction buffer (10X) "complete" with 25mM MgCl<sub>2</sub>  
MgCl<sub>2</sub> (100 mM)

Cat.-no	Description	Amount
N 140	Maximo DFS-Plus Taq Polymerase	500 units
N 142	Maximo DFS-Plus Taq Polymerase	5x500 units
N 144	Maximo DFS-Plus Taq Polymerase	20x500 units

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# Maximo M-SuperhotTaq DNA Polymerase (qPCR)

## Features:

Maximo M-Superhot Taq DNA Polymerase for qPCR is designed for Real-Time PCR and Hot-start PCR. A special inhibitor suppresses the reaction at room temperature until after the first denaturation step. This prevents primer-dimers and other artefacts. Using the enzyme there is no need to adjust the existing standard PCR protocol.

## Applications:

- Hot Start and real time PCR
- Multiplex PCR
- Amplification of complex genomic and cDNA templates
- no primer-dimers and other artefacts; inactive at room temperature
- short activation time for real time PCR
- enhanced PCR sensitivity

## Description:

Maximo M-Superhot Taq DNA for qPCR and Hot-Start-PCR is an optimized mixture of a high processive Taq DNA Polymerase and special inhibitors to Taq DNA for real time PCR. The enzyme is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-stranded specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94KD. It is developed for real time PCR or as basis enzyme for real time PCR diagnostics systems.

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## Features:

- Anti-Taq inhibits the reaction of Polymerases at room temperature
- Ideal as enhancer for PCR reactions

## Application:

Polymerase detergent/enhancer for:

- High specific amplification
- Multiplex amplification
- High sensitivity applications
- Low-copy number PCR

## Unit definition:

One unit is defined as the amount of antibodies required to blocks 50% activity of 1 µg of Polymerase

## Note:

The ratio units/mg of polymerase varies (up to factor of 10). The binding rate Anti-Taq depends strongly of the polymerization region (active epitopes) of the polymerase. The optimized mixture has to be found in empiric test and for every lot of polymerase, again.

**Concentration:** 2-10 mg/ml

## Storage Buffer:

20mM Tris-HCl (pH 7.0, at 22 °C);50 mM KCl;0.1mM EDTA; 50% glycerol

## Units/mg-ratio:

2300 units of the specific polymerase activity are equal to 1 mg of antibodies

Cat.-no	Description	Amount
S 131	Maximo m-Anti Taq	100 µg
S 132	Maximo m-Anti Taq	500 µg

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## (qPCR)

### Features:

p- Taq Polymerase for qPCR is designed for Real Time PCR. Polyclonal antibodies inhibit the reaction at room temperature until after the first denaturation step. This prevents primer-dimers and other artefacts. Using the enzyme there is no need to adjust the existing standard PCR protocol in your lab.

### Applications:

- Hot Start and real time PCR
- Multiplex PCR
- Amplification of complex genomic and cDNA templates
- no primer-dimers and other artefacts; inactive at room temperature
- short activation time
- enhanced PCR sensitivity

### Description:

Maximo p- Taq Polymerase is an optimized mixture of a high processive Taq Polymerase and polyclonal antibodies to Taq Polymerase. The enzyme is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-strand specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94 kDa.

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# Pfu/Psp DNA Polymerase

## Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors.

## Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis

## Description:

Pfu/Psp DNA polymerase, isolated from the archae bacteria Pyrococcus furiosus/species is a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** 5 u/μl

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nM of dNTPs into acid insoluble material in 30 minutes at 75°C.

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(ready-to-load)

## Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors. For visual control and for fast loading on the gel the enzyme contains a red dye and loading buffer for agarose electrophoresis.

## Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis
- PCR where visual control is needed

## Description:

Pfu/Psp DNA polymerase ready-to-load is isolated from the archae bacteria Pyrococcus f-species is a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu/Psp DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** 1 u/μl

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nM of dNTPs into acid insoluble material in 30 minutes at 75°C.

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(2X Pre-mix, ready to use)

## Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors. To reduce the risk of contamination, pipetting errors and to increase the repeatability of results the 2X-preMix contains an optimized mixture of enzyme, dNTP's and reaction buffer. Just add your template DNA and primers.

## Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis
- PCR where visual control is needed

## Description:

Pfu/Psp DNA polymerase 2X-preMix is isolated from the archae bacteria Pyrococcus f-species, a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu/Psp DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** Premix 2X (25µl per reaction)

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nM of dNTPs into acid insoluble material in 30 minutes at 75°C.

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# Maximo Tth Polymerase

## Features:

The thermostability and the reverse transcriptase (RT) activity of Tth DNA polymerase is useful in amplifying DNA from RNA templates that contain G-C-rich sequences or secondary structures since the elevated temperatures serve to denature the template RNA. Higher temperatures (in contrast to other enzymes for RT-PCR) also result in increased specificity of primer hybridization and extension. The concentration of RNA template for effective reverse transcription with Tth DNA polymerase should be higher if to compare with reverse transcription directed by Reverse Transcriptase (M-MuLV, AMV).

## Applications:

- PCR and RT-PCR
- cDNA synthesis

## Description:

Tth DNA Polymerase is a thermostable enzyme that replicates DNA at 74 °C and exhibits a half-life of 20 minutes at 95 °C isolated from eubacterium *Thermus thermophilus* strain HB8. Tth catalyzes the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of magnesium and the polymerization of nucleotides into DNA using an RNA template in the 5'→3' direction in the presence of manganese. The enzyme has a molecular weight of 94 000 daltons as estimated from the predicted amino acid sequence and exhibits 5'→3' exonuclease activity. Tth is recommended for use in PCR, RT-PCR, reverse transcription and primer extension reactions at elevated temperature.

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# UDG (Uracil-DNA Glycosylase)

## Applications:

- site-directed mutagenesis
- as a probe for protein-DNA interaction studies
- Glycosylase mediated single nucleotide polymorphism detection (GMPD)
- SNP genotyping
- Rapid and efficient cloning of PCR products
- Elimination carry-over contamination in PCR

## Description:

The Uracil-DNA Glycosylase (UDG, UNG) catalyzes the hydrolysis of the N-glycosylic bond between the uracil and sugar, leaving an apyrimidinic site in uracil-containing single or double-stranded DNA. Shows no activity on RNA. Molecular weight: 25.6 kDa monomer.

**Concentration:** 20 - 40 u/µl

## Unit definition:

One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [3H]-uracil in a 50 µl reaction containing 0,2 µg DNA (104 – 105 cpm/µg) in 30 minutes at 37°C.

## Storage Buffer:

20 mM Tris-HCl (pH 8.0); 50 mM NaCl; 1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol

## Reaction Buffer 10X

200 mM Tris-HCl (pH 8.0 at 25°C), 10 mM EDTA, 10 mM DTT. Incubation of 1X buffer at 37°C.

Cat.-no	Description	Amount
111-005	Uracil-DNA Glycosylase (UDG)	200 units
111-025	Uracil-DNA Glycosylase (UDG)	1000 units

# PCR-Mastermixes for Real-time PCR (qPCR)

## Colorless-Line

- SuperHot qPCR Master Mix RT..... 34-35
- SuperHot qPCR Master Mix HY ..... 36-37
- SuperHot qPCR Master Mix HS ..... 38-39

## Safe-Line

### with EvaGreen

- Eva-Master Mix E1 ..... 40-41
- Eva-Master Mix E2 ..... 42-43
- Eva-Master Mix E3 ..... 44-45
- Eva Master Mix E4..... 46-47

### without stain DLP (for dual label probes)

- Master Mix DLP1 ..... 48-49
- Master Mix DLP2 ..... 50-51
- Master Mix DLP3 ..... 52-53
- Master Mix DLP4 ..... 54-55

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# SuperHot qPCR Master Mix RT

## Features:

- The Master Mix (2X) offers both: high sensitivity and specificity
- Time saving ready-to-use qPCR Mastermix
- repeatable and reliable results
- efficient PCR for a wide range of template concentrations
- activation time < 3 min

## Applications:

- Realtime PCR and quantitative PCR e.g. with Sybr green, Eva green or probes
- High-throughput PCR
- Multiplex PCR
- Low copy targets PCR

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase (m-Superhot-Taq) for qPCR, dATP, dCTP, dGTP, dTTP, reaction buffer with stabilizers and enhancers, 1 Tube PCR-grade water, 1 Tube  $MgCl_2$

Cat.-no	Description	Amount
S240	qPCR Master mix RT (2x1,25 ml)	100 rcs / 50µl

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# SuperHot qPCR Master Mix HY

## Features:

- The Master Mix (2X) is optimized for "high yield" PCR results
- Time saving ready-to-use qPCR Mastermix
- repeatable and reliable results
- efficient PCR for a wide range of template concentrations
- activation time < 3 minutes

## Applications:

- Realtime PCR and quantitative PCR e.g. with Sybr green, Eva green or probes
- DNA Labeling with biotin or radioactive nucleotides
- Multiplex PCR
- Sequencing of double or single stranded DNA
- Low copy targets PCR

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity of the mix is achieved by an optimized hot-start polymerase and an optimized Buffer system. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dTTP, reaction buffer with stabilizers and enhancers for "High-Yield" results, 1 Tube PCR-grade water, 1 Tube MgCl<sub>2</sub>

Cat.-no	Description	Amount
S250	SuperHot qPCR Master mix HY (2x1,25 ml)	100 rcs / 50µl

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## Features:

- The Master Mix (2X) offers improved specificity
- Time saving ready-to-use qPCR Mastermix
- repeatable and reliable results
- efficient PCR for a wide range of template concentrations
- activation time < 3 minutes

## Applications:

- Realtime PCR and quantitative PCR
- DNA Labeling with biotin or radioactive nucleotides
- Multiplex PCR
- Sequencing of double or single stranded DNA
- Low copy targets PCR

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity of the mix is achieved by an optimized hot-start polymerase and an optimized Buffer system. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dTTP, reaction buffer with stabilizers and enhancers for "High-Specificity" results, 1 Tube PCR-grade water, 1 Tube MgCl<sub>2</sub>

Cat.-no	Description	Amount
S260	SuperHot qPCR Master mix HS (2x1,25 ml)	100 rcs / 50µl

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# Eva-Master Mix E1

## Features:

- The Master mix contains dUTP instead of dTTP
- It contains EvaGreen as fluorescent dye
- The qPCR / RT-PCR Mastermix E1 is ready-to-use and is optimized for high specificity and sensitivity
- because of optimized reaction buffer
- easy to us because ready-to-use Master Mix
- The Master Mix can be used with ROX as reference dye (1x concentrated)

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix contains dUTP instead of dTTP and allows an UNG (Uracil-N-Glycosylase) treatment at the onset of thermal cycling to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Master mix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, EvaGreen, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S150	qPCR / real-time PCR Mastermix	100 rcs / 50µl

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## Eva-Master Mix E2

### Features:

- The Master mix contains dUTP instead of dTTP
- The Master contains UNG (Uracil-N-Glycosylase)
- It contains EvaGreen as fluorescent dye
- The qPCR / RTD-PCR Master mix E2 is ready-to-use and is optimized for high specificity and sensitivity
- because of optimized reaction buffer
- easy to us because ready-to-use Master Mix
- The Master Mix can be used with ROX as reference dye (1x concentrated)

### Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

### Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions. It contains optimized UNG (Uracil-N-Glycosylase) for the treatment at the onset of thermal cycling.

**Concentration:** The Mastermix is 2x concentrated

### List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, EvaGreen, optimized reaction buffer with KCl and MgCl<sub>2</sub>, Uracil-N-Glycosylase, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S160	Real-time PCR Mastermix E2	100 rcs / 50µl

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## Features:

- The Master mix contains dUTP instead of dTTP
- The Mix contains ROX (500nM) as passive Reference dye (it provides a baseline in multiplex reactions)
- It contains EvaGreen as fluorescent dye
- The qPCR / RTD-PCR Master mix E3 is ready-to-use and is optimized for high specificity and sensitivity
- because of optimized reaction buffer
- easy to us because ready-to-use Master Mix

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, EvaGreen, ROX, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S170	Real time PCR Mastermix E3 (2x1,25ml)	100 rcs / 50µl

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## Features:

- The Master mix contains dUTP instead of dTTP
- The Master contains UNG (Uracil-N-Glycosylase)
- The Mix contains ROX (500nM) as passive Reference dye (it provides a baseline in multiplex reactions)
- It contains EvaGreen as fluorescent dye
- The qPCR / RTD-PCR Master mix E4 is ready-to-use and is optimized for high specificity and sensitivity
- because of optimized reaction buffer
- easy to us because ready-to-use Master Mix

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, EvaGreen, ROX, UNG, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S180	qPCR / Real-time PCR Mastermix E4	100 rcs / 50µl

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## Features:

- optimized realtime PCR Mastermix using probe based detection (e.g. FRET, Molecular Beacons or TaqMan)
- The Master mix contains dUTP instead of dTTP
- The qPCR / RT-PCR Mastermix DLP1 is ready-to-use and is optimized for high specificity and sensitivity because of optimized reaction buffer
- easy to use because ready-to-use Master Mix for block based PCR Cycler
- The Master Mix can be used with ROX as reference dye (1x concentrated)

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix contains dUTP instead of dTTP and allows an UNG (Uracil-N-Glycosylase) treatment at the onset of thermal cycling to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Master mix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S190	qPCR Master mix DLP1	100 rcs / 50µl

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## Features:

- optimized real time PCR Mastermix using probe based detection (e.g. FRET, Molecular Beacons or TaqMan)
- The Master mix contains dUTP instead of dTTP
- The Master contains UDG (Uracil-Glycosylase)
- The qPCR / RTD-PCR Master mix DLP is ready-to-use and is optimized for high specificity and sensitivity because of optimized reaction buffer
- easy to use because ready-to-use Master Mix
- The Master Mix can be used with ROX as reference dye (1x concentrated)

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions. It contains optimized UDG (Uracil-Glycosylase) for the treatment at the onset of thermal cycling.

**Concentration:** The Mastermix is 2x concentrated

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## Features:

- The Master mix contains dUTP instead of dTTP
- The Mix contains ROX (500nM) as passive Reference dye (it provides a baseline in multiplex reactions)
- The qPCR / RTD-PCR Master mix DLP3 is ready-to-use and is optimized for high specificity and sensitivity because of optimized reaction buffer
- easy to use because ready-to-use Master Mix

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Mastermix is 2x concentrated

## List of components qPCR / RTD-PCR Master mix:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, ROX, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S210	qPCR Master mix DLP3 (2x1,25ml)	100 rcs / 50µl

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## Features:

- The Master mix contains dUTP instead of dTTP
- The Master contains UDG (Uracil-Glycosylase)
- The Mix contains ROX (500nM) as passive Reference dye (it provides a baseline in Multiplex reactions)
- The qPCR / RTD-PCR Master mix DLP4 is ready-to-use and is optimized for high specificity and sensitivity because of optimized reaction buffer
- easy to us because ready-to-use Master Mix

## Applications:

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

## Description:

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup. The mix offer dUTP instead of dTTP and UDG to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Mastermix is 2x concentrated

## List of components:

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, ROX, UNG, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers, PCR-grade water

Cat.-no	Description	Amount
S220	qPCR Master mix DLP4 (2x1,25ml)	100 rcs / 50µl

## L-Mastermix for lyophilisation (low glycerol content)

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## L-Mastermix for lyophilisation (low glycerol content)

### Applications:

- Preparation of customized lyophilized PCR Master Mix
- Routine PCR
- Primer extension
- Low-copy target PCR
- Multiplex PCR
- Real-Time PCR (qPCR)

### Description:

5x L-Realtime PCR Master MIX is a ready-to-use premix of all components for "Real-time" amplification or standard PCR of target DNA. The cocktail contains stabilizers and enhancers, which improve thermo-stabilization of enzyme during PCR amplification and storage. The Mix is optimized for PCR with complex, low-copy number DNA templates, multiplex PCR, and improving specificity of PCR. 5x L-Realtime PCR Master MIX will be delivered in a customized concentration of  $MgCl_2$  and Hot-Start Enzyme content.

The mix is designed and prepared **for lyophilisation in tubes, strips, plates.**

### Stability:

5X L-Mix PCR Master MIX is stable for 24 months at -20°C, or for 6 months at +4°C storage without freezing in lyophilized form.

### Content: (1X) as an example:

- e.g. 2.5 Units m-SuperHot Taq DNA Polymerase (customized concentration)
- 0.2mM each of dNTP's
- e.g. 3 mM  $MgCl_2$  (customized concentration)
- Reaction Buffer: 16,6mM  $(NH_4)_2SO_4$ ; 67,0mM Tris-HCL (pH 8.8 at 25°C)
- Stabilizer/enhancer

**The final concentration of  $MgCl_2$  and Hot-Start Enzyme and in 5x stock-solution is variable according to the customers order**

Cat.-no	Description	Amount
S270	L-Master mix for qPCR	50 ml

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
# LightCycler Capillaries

## qPCR-Kapillaren

High precision polycarbonate capillaries are designed for use in LightCyclers from Roche. The polycarbonate capillaries are more stable in comparison with glass ones, thus the risk of breakage the capillaries is minimal. The cap provides a secure seal, minimizing the risk of cross contamination and reducing the risk of injuries, PCR protocols and handling procedures can be transferred. Without any changes from the original Roche manual in most experiments. PCR in PC capillaries shows sensitivity comparable to those in glass capillaries. The capillaries can hold reaction volumes from 10 to 50 µl. They are packed in a convenient 96 format, that ensures easy handling. The capillaries are shipped together with tide fitting caps that can be removed after PCR for recovery of PCR products.

## Transfer-Pin

For more comfortable handling of the capillaries and the caps a special transfer pin can be used. With this transfer pin you can take and hold the caps as well as the capillaries handle them as usual, but more comfortable.

# LightCycler Capillaries


## Carousel

This carousel for LightCyclers Version 1.2/1.5 or 2.0 from Roche is needed to enable the use of the capillaries. The carousel has 32 positions to hold the polycarbonate capillaries. With this carousel the capillaries are positioned precisely above the detection unit of the cyclor. The Rotor is manufactured ensuring maintenance of shape and dimensions during its product life. The Rotor is almost unbreakable. This carousel is part of our LightCycler Starter

### NOTE:

GeneOn assumes no liability for the improper use of the carousel in the centrifuge.

### Your advantages:

- no risk of breaking
- no risk of contamination
- same sensitivity
- no or less change in protocol
- about 40 to 50 % better price
- easy post sample recovery

Cat.-no	Description	Amount
A31-960	LightCycler qPCR-capillaries: in racks with Caps	10x96
A31-30SKO	LightCycler Starter-Set I : Carousel-Adapter for LightCycler 1,2 or 1,5 and 1 rack, transferpin and cap	1 Rotor + 1x96 + transferpin
A31-30SK1	LightCycler Starter-Set II: Carousel-Adapter for LightCycler 2 and 1 rack, transferpin and cap	1 Rotor + 1x96 + transferpin
A31-30TP	Transfer Pin	1 pcs

## PCR Mastermixes for standard PCR

- Maximo Taq DNA Polymerase (ready-to-use)..... 64-65
- Maximo Blue-Taq DNA Polymerase (ready-to-load) ..... 66-67
- Pfu RED DNA Polymerase (ready-to-load) ..... 68-69
- Pfu DNA Polymerase (2X Pre-mix, ready to use)..... 70-71

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(2X pre-mix, ready-to-use)

## Features:

Maximo Taq DNA Polymerase 2X-preMix provides robust PCR performance in a wide range of PCR applications and different templates. Best value in terms of cost per unit. The optimized mixture of all components reduces pipetting mistakes and ensures repeatable results - every day.

## Applications DFS-Taq Polymerase:

- Standard / General PCR
- optimized for high specificity
- High-throughput PCR, automated pipetting, or plate based PCR
- Gene mutation
- T/A cloning

## Description:

Maximo Taq DNA Polymerase 2X-preMix is optimized and ready-to-use mixture of all components for a successful PCR. Only your primers and your DNA Template has to be added. Maximo Taq DNA Polymerase 2X-preMix contains a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-strand specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94KD.

**Concentration:** the mixture is 2X concentrated

## Unit definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitation material in 30min at 74 degree

## List of components:

0.1U/ul Taq DNA Polymerase, 0.4 mM dATP, 0.4 mM dGTP, 0.4 mM dCTP, 0.4 mM dTTP, 4 mM MgSO<sub>4</sub>, 20 mM KCl, 16 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20 mM Tris-HCl, pH8.8

Cat.-no	Description	Amount
S 113	Maximo Taq DNA Polymerase 2X-preMix	1x100 rcs (2x1.25 ml)
S 114	Maximo Taq DNA Polymerase 2X-preMix	10x100 rcs (20x1.25 ml)

# Maximo Blue-Taq DNA Polymerase (ready-to-load)

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# Maximo Blue-Taq DNA Polymerase (ready-to-load)

## Features:

Maximo Taq-Blue DNA Polymerase provides robust PCR performance in a wide range of PCR applications and different templates. After PCR reaction the enzyme can be loaded to agarose gel, no dye and DNA loading buffer is needed. The enzyme is time- and cost saving, because it includes dye and loading buffer, already.

## Applications:

- Standard / General PCR with visible control
- High-throughput PCR
- Primer extension
- Gene mutation
- T/A cloning

## Description:

Maximo Tag-Blue DNA Polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double-strand specific 5'→3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94KD.

**Concentration:** 1 u/μl

## Unit definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitation material in 30min at 74 degree

## Storage Buffer:

25mM Tris-HCl (pH8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, 0.5% Nonident P 40, 0.5% Tween 20, in blue loading buffer

## Reaction Buffers supplied with the enzyme:

**10X Buffer I:** 500mM KCl, 100mM Tris-HCl, pH 9.0, 1% Triton X-100, 15mM MgCl<sub>2</sub>

**10X Buffer II:** 500mM KCl, 100mM Tris-HCl, pH 9.0, 1% Triton X-100

**MgCl<sub>2</sub>:** 100mM

Cat.-no	Description	Amount
S 111	Maximo Taq-Blue DNA Polymerase	500 units
S 112	Maximo Taq-Blue DNA Polymerase	5x500 units
S 128	Maximo Taq-Blue DNA Polymerase	20x500 units
S 129	Maximo Taq-Blue DNA Polymerase	100x500 units (or as Bulk)

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(ready-to-load)

## Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors. For visual control and for fast loading on the gel the enzyme contains a red dye and loading buffer for agarose electrophoresis.

## Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis
- PCR where visual control is needed

## Description:

Pfu/Psp DNA polymerase ready-to-load is isolated from the archae bacteria Pyrococcus f-species is a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu/Psp DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** (1 u/μl)

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nM of dNTPs into acid insoluble material in 30 minutes at 72°C.

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(2X Pre-mix, ready to use)

## Features:

Pfu/Psp DNA polymerase replicates DNA at 75°C catalyzing the polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>+</sup>. Pfu DNA polymerase possesses 3' to 5' exonuclease proof reading activity that enables the polymerase to correct nucleotide-misincorporation errors. To reduce the risk of contamination, pipetting errors and to increase the repeatability of results the 2X-preMix contains an optimized mixture of enzyme, dNTP's and reaction buffer. Just add your template DNA and primers.

## Applications:

- blunt end PCR cloning
- PCR and primer extension where "high fidelity" is required
- Site-directed mutagenesis
- PCR where visual control is needed

## Description:

Pfu/Psp DNA polymerase 2X-preMix is isolated from the archae bacteria Pyrococcus f-species, a thermostable Polymerase of approximately 90000 daltons. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. The Pfu/Psp DNA Polymerase has no detectable reverse transcriptase activity.

**Concentration:** Premix 2X (25µl per reaction)

## Unit definition:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nmol of dNTPs into acid insoluble material in 30 minutes at 75°C.

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# AMV Reverse Transcriptase

## Applications:

- RT PCR
- Synthesis of cDNA
- RNA Sequencing

## Description:

AMV Reverse Transcriptase (AMV RT) catalyzes the polymerization of DNA using template DNA, RNA or RNA:DNA hybrids. The enzyme possesses an intrinsic RNase H activity. AMV RT possesses multiple enzymatic activities including RNA- and DNA-directed DNA polymerase, DNA-RNA unwinding activity, a sequence-specific  $Mn^{2+}$ -dependent endonuclease and ribonuclease H.

**Concentration:** 10 u/ $\mu$ l

## Unit definition:

One unit is the amount of enzyme required to catalyze the transfer of 1nmol of deoxynucleotide into acid-precipitable material in 10 minutes at 37°C.

## Storage Buffer:

200mM potassium phosphate (pH7.2), 0.2% Triton X-100, 2mM DTT and 50% glycerol

## Reaction Buffer 5X:

250mM Tris-HCl (pH8.3 at 25°C), 250mM KCl, 50mM  $MgCl_2$ , 2.5mM spermidine and 50mM DTT

Cat.-no	Description	Amount
105-400	AMV Reverse Transcriptase	200 units
105-410	AMV Reverse Transcriptase	5x200 units

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## Applications:

- Inhibits common eukaryotic RNases
- Active over a broad pH range (pH 5-8)
- high levels of inhibition over a wide range of conditions

## Description:

RNase Inhibitor is a recombinant human placental protein which inhibits ribonucleases (RNases) A, B and C. It does not inhibit RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from Aspergillus. There is no inhibition of polymerase activity when the protein is used with Taq DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

**Concentration:** 40 u/µl

## Unit definition:

One unit is the amount of enzyme required to inhibit by 50% the activity of 5 ng of RNase A at 25°C (This inhibitor activity is determined by its ability to inhibit hydrolysis of cyclic 2', 3'-CMP by RNase A).

## Storage Buffer:

20 mM HEPES-KOH (pH7.6), 50 mM KCl, 8 mM DTT and 50% glycerol

Cat.-no	Description	Amount
105-310	Ribonuklease (RNase) Inhibitor	2000 units
105-350	Ribonuklease (RNase) Inhibitor	10000 units

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## Reverse (M-MLV RT)

### Applications:

- RT PCR
- Synthesis of cDNA
- mRNA 5'-end Mapping by Primer Extension Analysis
- End-labeling of DNA
- Dideoxynucleotide Sequencing

### Description:

MMLV Reverse Transcriptase, encoded by Moloney Murine Leukemia Virus (MMLV RT) is an RNA-dependent DNA polymerase that synthesizes the cDNA first strand from a single-stranded RNA template to which a primer has been hybridized. MMLV RT will also extend primers hybridized to single-stranded DNA. Second strand cDNA synthesis can be achieved from some mRNA templates without an additional DNA polymerase.

**Concentration:** 200 u/µl

### Unit definition:

One unit of the enzyme incorporates 1 nmol dTTP into acid-precipitable material in 10 minutes at 37°C, using poly(A) oligo dT as a template primer.

### Storage Buffer:

200 mM potassium phosphate (pH 7.2), 0.2% Triton X-100, 2 mM DTT and 50% glycerol

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# Maximo Tth Polymerase

## Features:

The thermostability and the reverse transcriptase (RT) activity of Tth DNA polymerase is useful in amplifying DNA from RNA templates that contain G-C-rich sequences or secondary structures since the elevated temperatures serve to denature the template RNA. Higher temperatures (in contrast to other enzymes for RT-PCR) also result in increased specificity of primer hybridization and extension. The concentration of RNA template for effective reverse transcription with Tth DNA polymerase should be higher if to compare with reverse transcription directed by Reverse Transcriptases (M-MuLV, AMV).

## Applications:

- PCR and RT-PCR
- cDNA synthesis

## Description:

Tth DNA Polymerase is a thermostable enzyme that replicates DNA at 74 °C and exhibits a half-life of 20 minutes at 95 °C isolated from eubacterium *Thermus thermophilus* strain HB8. Tth catalyzes the polymerization of nucleotides into duplex DNA in the 5' -->3' direction in the presence of magnesium and the polymerization of nucleotides into DNA using an RNA template in the 5' -->3' direction in the presence of manganese. The enzyme has a molecular weight of 94 000 daltons as estimated from the predicted amino acid sequence and exhibits 5' -->3' exonuclease activity. Tth is recommended for use in PCR, RT-PCR, reverse transcription and primer extension reactions at elevated temperature.

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## Oligo(dT)<sub>15</sub> (1DA)

### Description:

- The primer is designed to initiate synthesis of a cDNA from total RNA in a reverse transcription reaction.
- For use in generating labeled cDNA for screening microarrays.
- The primer hybridizes to the poly(A) tail of mRNA.

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## Random Primers

### Description:

The primer is a mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Random hexamer Primers can be used for first-strand cDNA synthesis and cloning. The primer is supplied in lyophilized form.


Cat.-no	Description	Amount
S 300	Random Primer	30 µg

## Nucleotides Ultrapure

- Nucleotide-Mix 40mM (4x10mM each) /  
Nucleotide Set 100mM (4 x 1 ml) ..... 88-89
- Biotin-11-dUTP (1mM)..... 90-91

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
## Applications:

- all molecular biology applications
- PCR/qPCR
- reverse transcription
- DNA labeling
- DNA sequencing

## Description:

dNTP-sets from GeneON have PCR Grade and contains four separate tubes of dATP, dCTP, dGTP and dTTP supplied as aqueous solutions at pH 8.5. The dNTP-Mix contains an optimized mixture of dNTP's with 10 mM of each nucleotide.

## Quality control:

- purity: purity test with RP-HPLC: greater 99 % (stable from lot-to-lot)
- Amplification tests:  
 Lambda DNA: 18 kb fragment, less than 20 pg template  
 RT-PCR: Human DNA, 600 bp fragment, less than 20 pg template
- Free of bacterial - and human DNA
- Free of: DNase, RNase, Protease and no nicking activity
- produced in German factory (ISO: 9001/2001)

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## Biotin-11-dUTP (1mM)

### Applications:

- labeling of DNA
- nick translation of DNA
- 3'-end labeling
- cDNA synthesis or primer extension reactions
- incorporation into DNA using: Reverse Transcriptase, Taq DNA Polymerase, Klenow Fragment

### Description:

Biotin-11-dUTP (Biotin-11-2'-deoxyuridine-5'-triphosphate, tetralithium salt) compound for non-radioactive DNA labeling.

The number "11" is the number of carbon atoms in the backbone of linker between dUTP and biotin. The longer the linker is, the more effective interaction of biotin with avidin occurs. The length of spacer "11" is optimal for most of applications.

**Quality:** more than 96% (IEC, TLC, NMR, UV)

**Concentration:** 1 mM

### Storage buffer:

10 mM TrisHCl, pH 7.5, 1 mM EDTA

Cat.-no	Description	Amount
110	Biotin-11-dUTP	100 µl
111	Biotin-11-dUTP	5x100 µl

● GeneON ●